



## Automated DNA Extraction and Species Analysis from Cheese

### Introduction

People are more and more aware of the fact that nutrition has a significant impact on the overall quality of life. Additionally, food scandals of the recent past even sharpened this awareness. For those reasons, customers expect stringent and reliable controls of the ingredients used for food production. Some of these controls (e.g. species identification, allergen identification, halal testing) can be performed using molecular biological methods. The general workflow for molecular biological analysis consists of nucleic acid extraction from sample material followed by quality analysis (e.g. by photometry) and a diagnostic assay mostly based on polymerase chain reaction (PCR). The extraction of nucleic acids from food is challenging due to extensive processing and usage of additives. As soon as nucleic acids are extracted the analysis is dependent on sufficient quantity and quality of DNA and the availability of an assay to detect the parameter of interest, e.g. for species identification.

Herein we describe extraction of deoxyribonucleic acid (DNA) using InnuPure® C16 touch with innuPREP Food DNA Kit-IPC16 specifically developed for extraction of DNA from processed food. After quality control using ScanDrop® 250 and CHIPCUVETTE®, extracted DNA is analyzed for the presence of sheep-, goat- and cow-specific genes using innuDETECT Cheese Assay and qTOWER<sup>3</sup> G (qPCR) system. The workflow for food quality control presented here allows automated extraction of up to 16 food samples in parallel. This enables medium sample throughput with a minimum hands-on time and manual operations. Downstream analysis includes parallel detection of DNA quantity and quality with CHIPCUVETTE® microfluidic system and species analysis with highly specific and sensitive qPCR.

### Challenge

Analysis of the source species of milk for cheese production.

### Solution

Combination of automated nucleic acid extraction, photometry and quantitative real-time PCR (qPCR) for easy-to-use and highly sensitive species analysis for cheese.

## Materials and Methods

### Samples and Reagents

- Two types of commercially available cheese, either declared as derived from goat milk or cow milk
- innuPREP Food DNA Kit-IPC16
- innuDETECT Cheese Assay
- 96 Well PCR-Plate (0.2 ml; LP), full-skirted, white
- Optical sealing foil (77 × 140 mm), adhesive, transparent, peelable

### Instrumentation

- InnuPure® C16 *touch*
- BioShake iQ with adapter for 35 × 1.5 ml tubes or 24 × 0.5 ml tubes
- ScanDrop® 250 with CHIPCUVETTE®
- qTOWER<sup>3</sup> G

### Preparation of DNA from cheese

DNA is extracted with innuPREP Food DNA Kit-IPC16 and InnuPure® C16 *touch* according to the following protocol.

200 mg of the cheese sample is transferred into 1.5 ml reaction tubes. The sample is mixed with 800 µl Lysis Solution CBV and 20 µl proteinase K, followed by vigorous mixing on a vortex. For homogenization, lysis and proteinase K digest the sample is incubated for 1 h at 65 °C shaking with 1000 rpm. Subsequently, solid particles are separated by centrifugation of the lysate for 10 min with 11000 × g. 400 µl of the supernatant is transferred into the Reagent Plate of the innuPREP Food DNA Kit-IPC16. The sample tray carrying the Reagent Plate with samples, disposable pipet tips and elution tubes is put into InnuPure® C16 *touch*, followed by start of the protocol (Internal Lysis 200 µl – 04) of the InnuPure® C16 *touch*. The end of the extraction protocol has to be confirmed and the extracted DNA can be unloaded. DNA is stored at -20° C until further use. Each cheese variety is extracted in quadruplicates

### Analysis of DNA quantity and quality

The DNA is quantified using ScanDrop® 250 and CHIPCUVETTE®. The ScanDrop® family of instruments are polychromatic photometers assessing a whole spectrum within a single and quick detection process. The software for device control is preset to determine the quantity of DNA based on the absorption value of light with a wavelength of 260 nm. The quality is determined by calculation of the ratios of the absorption of light with wavelengths of 260 nm and 280 nm.



Figure. 1: CHIPCUVETTE®

Analytik Jena's CHIPCUVETTE® (Figure. 1) allows analysis of up to 16 samples in parallel (1 reference and 15 samples). Each sample can be detected with pathlengths of 1 mm and 0.1 mm making dilution redundant.

### Quantitative real-time PCR analysis of cheese samples

The innuDETECT Cheese Assay contains three individual primer/probe mixes for detection of beef-, sheep and goat-specific genes. Extracted DNA from each sample was mixed with each of the primer/probe mixes according to table 1.

Table 1: Pipetting scheme for innuDETECT Cheese Assay.

Component	V [ $\mu$ l]
2 $\times$ Mastermix	10
primer/probe mix	3
internal control	1
sample	1
PCR-grade H <sub>2</sub> O	ad 20

For non-template control (NTC), PCR-grade water is used instead of sample.

96 well plates and optical sealing films were used as described above. The qTOWER<sup>3</sup> G was programmed as follows (table 2).

Table 2: Temperature profile for qPCR.

T [ $^{\circ}$ C]	t [s]	Cycles
95	120	
95	10	35
62	45	

The device control software qPCRsoft is used for analysis and data presentation.

### Results and Discussion

After extraction of DNA from 200 mg of cheese, concentration and quality were analyzed using photometry (table 3).

Table 3: Results of photometric analysis of extracted DNA using ScanDrop<sup>®</sup> 250.

Cheese type	Replicate	$A_{260}/A_{280}$	c [ng/ $\mu$ l]
Goat	1	2.03	9.86
	2	2.02	11.07
	3	2.01	12.56
	4	2.01	12.59
Cow	1	2.07	16.3
	2	1.96	12.78
	3	2.03	12.38
	4	2.01	9.84

The photometric data of table 3 show that 200 mg cheese is sufficient for the cheese varieties tested to yield sufficient DNA for PCR-based downstream applications. Moreover, purity of DNA is good, although a complex and highly processed sample material was used.




Figure 2: Detection of sheep-specific target gene



Figure 3: Detection of goat-specific target gene.




Figure 4: Detection of cow-specific target gene.

qPCR analysis of the samples using the primer/probe mix for sheep-specific genes is negative for all samples of the two cheese varieties (figure 2). This correlates with expectations as the cheese is declared to be derived from cow and goat milk, respectively. According to this, the quadruplicates of the cheese declared as goat cheese are positive for detection of the goat-specific genes and negative for cow-specific genes (figure 3) while the cheese declared to be derived from cow milk is negative for goat-specific genes and positive for cow-specific genes (figure 4). False negative results can be excluded as Internal Positive Control is detected in DNA of all extracted samples.

## Conclusion

The results show that the workflow of DNA extraction using InnuPure® C16 *touch* with innuPREP Food DNA Kit-IPC16 in combination with innuDETECT Cheese Assay and qTOWER<sup>3</sup> allows safe and reliable determination of the origin of milk in cheeses.