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Contents

1	Basic Information	25
1.1	User manual notes	25
	User manual conventions	25
1.2	Intended use	25
1.3	Warranty and liability	26
2	Technical data	27
3	Safety instructions	28
3.1	General notes	28
3.2	Symbols and signal words used in the user manual	28
3.3	Safety marks of the UV lamp	29
3.4	Technical condition	30
3.5	Requirements for the operating personnel	30
3.6	Safety instructions for transport and installation	30
3.7	Safety instructions for operation	31
3.8	Handling hazardous substances	31
3.9	Chemical resistance of the device	31
3.10	Behavior during emergencies	32
3.11	Standards and directives	32
4	Technical description	33
4.1	Fundamental layout and function	33
4.2	Design, connections, indications	33
5	Site requirements	36
5.1	Installation conditions	36
5.2	Space requirement	36
5.3	Energy supply	36
6	Startup/operation	37
6.1	Installing the UV-lamp and switching it on	37
6.2	Switching off the UV lamp	38
6.3	Setting the operating time of the UV lamp with the time relay	39
	Operating time in dependence of coding & setting value	39
7	Maintenance and care	40
7.1	Cleaning the UV lamp	40
7.2	Changing the fuse	41
8	Transport and storage	42
8.1	Transport	42
8.2	Storage	42
9	Disposal	43

Figures

Fig. 1	Safety marks on front and top of the UV lamp.....	29
Fig. 2	UV-C radiation killing microorganisms [1].....	33
Fig. 3	Total view of the UV lamp for the InnuPure C16 <i>touch</i>	34
Fig. 4	Time relay on the right side of the device.....	34
Fig. 5	Safety switch on the left side of the device.....	34
Fig. 6	Inserting lamp into InnuPure C16 <i>touch</i> , connecting and switch-on	38
Fig. 7	Time relay on the right side of the device.....	39
Fig. 8	Operating time of the UV lamp at different settings.....	39

1 Basic Information

1.1 User manual notes

The InnuPure C16 *touch* UV lamp is intended for operation by qualified specialist personnel observing this user manual.

The user manual informs about the design and operation of the UV lamp for the InnuPure C16 *touch* and provides personnel familiar with the isolation and purification of nucleic acid the necessary know-how for the safe handling of the UV lamp. The user manual further includes notes on the maintenance and service of the UV lamp and potential causes and remedies of any faults.

User manual conventions

Instructions for action which occur in chronological order are numbered and combined into action units and furnished with the corresponding results.

Lists which are not in chronological order are shown as itemized lists, sub-listings as bullet points.

Safety notes are indicated by pictographs and signal words. The type and source of the danger are stated together with notes on preventing the danger. The meaning of the pictographs and signal words used is explained in the chapter "Safety instructions".

1.2 Intended use

The UV lamp is an optional accessory for the InnuPure C16 *touch* and is used for the internal contamination. The InnuPure C16 *touch* is a flexible and efficient extraction system for the fully automated isolation and purification of nucleic acid in chemical and biological laboratories. The equipment is intended for general laboratory use.

The UV lamp for the InnuPure C16 *touch* must only be used for the applications described in this instruction manual. The manufacturer does not accept liability for any other use, including that of any individual modules or components.

- The UV lamp must only be operated by qualified and trained personnel.
- The UV lamp must only be used in accordance with this instruction manual. This applies in particular to the adherence to the connection values, conditions of use and notes on maintenance, service, transport and disposal.
- The safety notices in this instruction manual must be observed.

It is not permissible

- to use the UV lamp outside the InnuPure C16 *touch*
- to use the UV lamp when the upper front door of the InnuPure C16 *touch* is open.

The operator is responsible for the use of the equipment as intended.

1.3 Warranty and liability

The warranty and liability period comply with the legal requirements and the provisions in the general terms and conditions of Analytik Jena AG.

Deviations from the intended use described in this user manual result in limitations of warranty and liability during a damage event. Damage to wearing parts is not included in the warranty.

Warranty and liability claims are excluded for personal injury and property damage due to one or several of the following causes:

- usage of the UV lamp for the InnuPure C16 *touch* not according to the intended purpose
- improper commissioning, operation and service of the UV lamp
- modifications of the UV lamp without prior consultation with Analytik Jena AG
- unauthorized handling inside the UV lamp
- operation of the UV lamp with faulty safety equipment or improperly fitted safety and protection equipment
- inadequate monitoring of the equipment components subject to wear
- use of other than original spare parts or wearing parts
- improper repairs
- faults due to the non-observance of this user manual

2 Technical data

System parameters	
Decontamination principle	UV radiation
Decontamination duration	Selectable from 0.05 s to 100 h
	Factory setting 30 min
Maintenance	maintenance-free
Cleaning of the UV lamp	Little contact surfaces with the interior of the InnuPure C16 <i>touch</i> . The frame is easy to disinfect.
Other technical data	
Weight	Approx. 1.5 kg
Dimensions (W x H x D)	263 mm x 129 mm x 450 mm
Electrical variables	
Power supply	230 V AC, 50 Hz, optionally 110 – 115 V AC, 60 Hz
Protection	1 x 1 AT / 250V
Max. power consumption	30 W
Environmental conditions	
Temperature during operation	15 – 35 °C
Humidity during operation	Max. 80 % at +30 °C
Temperature during storage	-10 – 55 °C
Humidity during storage	Max. 80 % (use desiccant)
Protection type	IP 20

3 Safety instructions

3.1 General notes

For your own safety and to ensure an error-free and safe operation of the UV lamp for the InnuPure C16 *touch*, please read this chapter carefully before using the appliance.

Follow all safety instructions described in this user manual.

You must also comply with the safety instructions and the information concerning handling, storage and disposal of the InnuPure C16 *touch*.

Besides the safety instructions in this user manual and the local safety regulations that apply to the operation of the device, the general applicable regulations regarding accident prevention, occupational health and safety and environmental protection have to be observed and complied with.

References to potential dangers do not replace the work protection regulations which must be observed.

3.2 Symbols and signal words used in the user manual

The user manual uses the following symbols and signal words to indicate hazards or instructions. The safety instructions are always placed before an action.

WARNING

Indicates a potentially hazardous situation which might result in death or very serious injuries (deformities) if not prevented.

CAUTION

Indicates a potentially hazardous situation which might result in light or minor injuries and material damage if not prevented.

ATTENTION

Provides information on potential material or environmental damage.

3.3 Safety marks of the UV lamp

Safety and notice symbols have been attached to the UV lamp for the InnuPure C16 *touch* whose content must always be observed.

Damaged or missing safety and notice symbols can cause incorrect actions leading to personal injury or material damage! The symbol labels must not be removed! Damaged symbol labels must be replaced without delay!

Fig. 1 Safety marks on front and top of the UV lamp

Warning symbol	Meaning
	Warning against a hazard location
	Warning of UV radiation
<p>WARNUNG: UV-Strahlen sind gefährlich für die Augen und die Haut! Benutzung nur im InnuPure C16 bei geschlossener Tür erlaubt!</p> <p>WARNING: UV-rays are dangerous for eyes and skin! Use only in the InnuPure C16 with closed door!</p>	Warning of UV radiation, intended use only inside the InnuPure C16 <i>touch</i> with closed door.
Notice symbol	Meaning
	Observe the instructions for use. Prior to starting work read the instructions for use.

3.4 Technical condition

The UV lamp for the InnuPure C16 *touch* corresponds in its design and construction to the current state of the art technology. Unauthorized modifications or changes, especially such that affect the safety of the staff and the environment, are generally not allowed.

- The operator must only operate the UV lamp in a sound and operationally safe condition. The technical condition must always comply with the legal requirements and regulations.
- Prior to every use, the UV lamp must be checked for damage and sound condition.
- Any changes in the UV lamp affecting its safety must be reported by the operating personnel to the operator without delay.

3.5 Requirements for the operating personnel

- The UV lamp for the InnuPure C16 *touch* must only be commissioned, operated and serviced by trained personnel instructed in technical safety.
- The operation or servicing of the UV lamp by minors or individuals under the influence of alcohol, drugs or medication is not permitted.
- It must be ensured that only authorized personnel works at the UV lamp.
- The operating personnel must be familiar with the dangers arising from the substances used inside the InnuPure C16 *touch*. The appropriate protective equipment must be used.

3.6 Safety instructions for transport and installation

- Only transport the UV lamp in its original packaging! Ensure that the UV lamp has cooled down and is cleaned and decontaminated before transport.
- To prevent health damage the following must be observed when moving the UV lamp in the laboratory (lifting and carrying):
 - Prevention of bumping or falling down since this could result in damage of the device or destruction of the illuminant. Beware of glass breakage.

3.7 Safety instructions for operation

- The operator of the UV lamp for the InnuPure C16 *touch* must make sure that the condition of the device including the safety equipment is sound before each commissioning. This applies in particular after each modification or extension of the UV lamp or its repair.
- The UV lamp must only be operated if all protective equipment (fuse) is in place, properly installed and fully operational.
- Free access to the mains plug in the power outlet must be ensured during operation.
- The UV lamp for the InnuPure C16 *touch* must not be removed from the device in switched-on state. The illumination area of the lamp must not be touched after switching off the lamp. Otherwise, the user might be burned (max. 60 °C).
- The mains plug must only be plugged into a power outlet with protective ground connection. The UV lamp may only be connected to power sources whose nominal voltage is the same as that on the type plate of the equipment. The protective effect must not be invalidated by the use of an extension line which does not have a protective conductor.
- Do not introduce any objects into the UV lamp and prevent fluids from entering the UV lamp through openings or gaps.
- Do not short-circuit the UV lamp and only use the fuse in accordance with the specifications in the technical data.

3.8 Handling hazardous substances

Even with intended use there is a risk of health damage when handling hazardous substances. The operator is solely responsible for the compliance with all safety requirements to protect individuals and property when handling radioactive, infectious, toxic, caustic, flammable and other hazardous substances.

- Control the handling of hazardous substances in accordance with the safety category of the lab, the details in the manufacturer safety data sheets and additional national and international regulations (WHO, "Laboratory Biosafety Manual").
- Wear personal protective equipment when working at the UV lamp.
- Observe all notices on the cleaning and decontamination of the UV lamp.

3.9 Chemical resistance of the device

Aggressive substances may damage the device. Although the materials used are resistant to most of the commonly used substances, material damage from aggressive substances may not be precluded.

- Before using any aggressive substances (e.g. acids, bases or organic solutions), check whether the materials in direct contact with the substances are resistant to them.
- Only use substances compatible with the materials listed.
- If in doubt contact Analytik Jena AG.

The housing parts of the UV lamp are not resistant to the following substances:

- | | |
|---|---|
| <input type="checkbox"/> hydrogen fluoride (HF/hydrofluoric acid) | <input type="checkbox"/> naphtha |
| <input type="checkbox"/> highly concentrated acids | <input type="checkbox"/> gasoline |
| <input type="checkbox"/> cleaning powder | <input type="checkbox"/> acetone |
| <input type="checkbox"/> paint thinner | <input type="checkbox"/> cleaning spray |

3.10 Behavior during emergencies

- Immediately unplug the mains plug from the mains outlet in situations of danger or in the event of accidents!

Because a rapid response can save lives during an emergency, the following has to be ensured:

- The operating staff must be familiar with the location of safety equipment, accident and danger alarms as well as first aid and rescue equipment as well as their handling.
- The operator is responsible for the training of the operating staff.
- All equipment for first aid (first-aid kit, eyewash bottles, stretcher, etc.) as well as equipment for firefighting (fire extinguishers) must be within reach and easy to access. All equipment has to be in a sound condition and has to be checked regularly.

3.11 Standards and directives

The UV lamp was manufactured according to the currently valid technology regulations and the approved safety related regulations.

When constructing the UV lamp, the relevant safety and health requirements of the applicable laws, standards and regulations were applied. The safety of the device is confirmed by the CE mark and the declaration of conformity.

Information regarding safety corresponds to the currently valid regulations of the European Union. In other countries the applicable laws and country specific regulations have to be complied with.

Besides the safety instructions in this user manual and the local safety regulations that apply to the operation of the device the general applicable regulations regarding accident prevention, occupational health and safety and environmental protection have to be observed and complied with.

References to potential dangers do not replace the work protection regulations which must be observed.

4 Technical description

4.1 Fundamental layout and function

The UV lamp is an optional accessory for an easy decontamination of the interior of the InnuPure C16 *touch* extraction systems after a purification of nucleic acids by means of UV radiation. The InnuPure C16 *touch* is a flexible and efficient extraction system for the fully automated isolation and purification of nucleic acid.

Working with molecular biological technologies often requires a contamination-free workplace. UV sterilization lamps provide protection against possible contamination in research laboratories, food-production enterprises, or in medical environments.

The installed UV lamp emits a maximum of UV radiation at a wavelength of 254 nm what is very effective for the destruction of bacteria, fungi, yeasts, and viruses. In addition, the UV lamp denatures unwanted DNA and RNA inside the InnuPure C16 *touch*.


Fig. 2 UV-C radiation killing microorganisms [1]

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4.2 Design, connections, indications

Beside the mains connection, the device front provides the button to switch on the UV lamp. During operation of the UV lamp, the green status LED is lit. The time relay for setting the operation time is located on the right side of the device. The factory setting is 30 minutes. A safety switch for monitoring the device environment is located on the left side of the device. The switch monitors the position of the UV lamp inside the InnuPure C16 *touch* and unlocks the time relay at contact of the particular interior wall of the device. With its frame the UV lamp can be positioned in the drawer for the sample tray. The UV lamp has little contact surfaces with the interior of the InnuPure C16 *touch*. The frame is easy to disinfect.

The type plate is attached to the bottom of the device. The type plate shows among other information the serial number of the device and the electrical connection data:

Information on the type plate	
Manufacturer (with address)	Analytik Jena, Konrad Zuse Str. 1, D-07745 Jena, Made in USA
Device type	UV Lamp InnuPure C16 touch
Order number	REF 37-400-096
Serial number	SN XXXX
Energy supply	230 V AC, max. 30 W / 50 Hz, 1 x 1 AT / 250 V
	110 – 115 V AC, max. 30 W / 60 Hz, 1 x 1 AT / 250 V
CE marking	
China RoHS label	The device contains substances. Analytik Jena AG guarantees that the substances will not leak within the next 25 years and, therefore, will not cause any environmental danger or health risk within this period of time if the device is used as intended.
Waste disposal symbol acc. to WEEE directive	Meaning: Do not dispose of as domestic waste!
Protection class of the housing	IP 20
Safety symbols	Meaning: Attention! Observe accompanying documents!
Year of manufacture	e.g. 2017
Device number	GN ZZZZ

5 Site requirements

5.1 Installation conditions

The following requirements apply to the climatic conditions in the operating room of the UV lamp (corresponding to those of the InnuPure C16 *touch*):

- Temperature range: +15 °C to +35 °C
- Max. humidity: 80 % at 30 °C

5.2 Space requirement

The space requirement results from the device dimensions:

W x H x D	263 mm x 129 mm x 450 mm
-----------	--------------------------

5.3 Energy supply

WARNING!

The UV lamp for the InnuPure C16 *touch* must only be connected to a properly grounded mains outlet in accordance with the voltage specifications on the type plate!

The UV lamp for the InnuPure C16 *touch* is operated on single-phase alternating current. The installation of the electrical equipment of the laboratory must comply with the standard DIN VDE 0100. After the connection point, an electrical current in accordance with the standard EN 60038 must be available.

6 Startup/operation

CAUTION!

Do not manually open the lower front door of the InnuPure C16 *touch*. Always open and close the lower front door via IP*extract* software. Do not reach inside the InnuPure C16 *touch* during operation. There is a risk of crushing and burning.

The InnuPure C16 *touch* must be switched off before the UV lamp is inserted into the device. During operation of the UV lamp the upper front door of the InnuPure C16 *touch* must be closed.

The UV lamp must be completely slid into the InnuPure C16 *touch* before the mains plug may be connected to the mains socket.

CAUTION! Danger of short-circuit from condensate!

If the storage and the installation temperature differ a lot, wait until the UV lamp has adapted to the new ambient temperature before connecting it in order to prevent damages to the device by condensation water.

ATTENTION

Retain the transport packaging! Return transport for service must be in the original packaging. This alone prevents transport damage.

6.1 Installing the UV-lamp and switching it on

1. Switch on the InnuPure C16 *touch* and the tablet.
The IP*extract* software starts automatically.
2. Select the item OPTIONS ▶ TOOLS ▶ DECONTAMINATION. Follow the prompts.
The IP*extract* software initializes the InnuPure C16 *touch*. The lower front door is automatically opened and the drawer with the sample is ejected.
3. Remove the sample tray from the InnuPure C16 *touch*.
4. Set the desired operating time via the time relay (see chapter "Setting the operating time of the UV lamp with the time relay" p. 39).
5. Position the UV lamp with its frame in the drawer for the sample tray.
6. Shut down the tablet via IP*extract*. Switch off the InnuPure C16 *touch* using the mains switch on the back of the device.
7. Completely slide the UV lamp into the InnuPure C16 *touch*.
8. Connect the mains plug to the UV lamp and to the mains socket.
9. Switch on the UV lamp via the front button and check whether the UV lamp gives light.
 - ✓ **The status LED is lit green during the set operating time.**

Fig. 6 Inserting lamp into InnuPure C16 *touch*, connecting and switch-on

6.2 Switching off the UV lamp

1. Wait until the set operating time has elapsed and the UV lamp is switched off.
 - ✓ **The green status LED at the front of the UV lamp turns off when the set operating time of the time relay has elapsed.**
2. Disconnect the mains plug from the UV lamp and the mains socket.
3. Slide the drawer with the UV lamp out of the InnuPure C16 *touch*. Remove the UV lamp.
4. Clean and disinfect the contact areas between the UV lamp and the interior of the InnuPure C16 *touch*.
5. Put the sample tray back into the InnuPure C16 *touch*.
6. Switch on the InnuPure C16 *touch* and the tablet.
7. Select the item **OPTIONS ▶ TOOLS ▶ INITIALIZATION**.

The drawer with the sample tray is retracted. The lower front door closes automatically.

CAUTION! Warming of the housing in the illumination area!

Depending on the operating time, the housing of the UV lamp may be heated up to 60° in the illumination area.

6.3 Setting the operating time of the UV lamp with the time relay

Depending on the coding of the DIP switches (1, 2, 3), the operating time of the UV lamp can be set with the dial. The time can be set from 0.05 s to 100 h. As an example, the table displays all possible coding for four setting values of the time relay. The factory setting of the operating time is 30 min.


Fig. 7 Time relay on the right side of the device

Operating time in dependence of coding & setting value

Setting value of: Time relay	DIP switch coding						
	(0.05...1)s	(0.5...10)s	(5...100)s	(0.5...10)min	(5...100)min	(0.5...10)h	(5...100)h
0.5 (minimum)	0.05s	0.,5s	5s	0.5min	5min	0.5h	5h
1	0.1s	1s	10s	1min	10min	1h	10h
5	0.5s	5s	50s	5min	50min	5h	50h
10 (maximum)	1s	10s	100s	10min	100min	10h	100h

Fig. 8 Operating time of the UV lamp at different settings

7.2 Changing the fuse

WARNING! Electric shock!

Disconnect the mains plug before changing fuses!

If the fuse of the UV lamp is faulty, you can replace it. The fuse holder comprises an additional fuse for the first exchange.

1. Disconnect the mains plug from the mains socket.
2. Open the fuse holder by pulling its cover (5 in Fig. 3 p. 34).
3. Replace the faulty mains fuse (rear fuse). Use the following fuse:
1 x 1 AT / 250 V
4. Close the fuse holder.

8 Transport and storage

8.1 Transport

ATTENTION!

Only use the transport packaging recommended by Analytik Jena AG for packaging. Only this provides optimum protection of the UV lamp against transport damage.

Observe the safety instructions in section "Safety instructions for transport and installation" p. 30. Transport the UV lamp very carefully to prevent damage from impact or vibration. The device should be transported in such a way that major temperature fluctuations are avoided and the formation of condensate is thus prevented.

8.2 Storage

ATTENTION

Environmental influences and condensate formation can destroy individual components of the UV lamp! The storage atmosphere of the UV lamp must be low in dust and free from aggressive vapors.

If the device is not installed immediately after delivery or is not required for prolonged periods, it should be stored in its original packaging. A suitable desiccant should be added to the equipment to prevent damage from moisture.

The following storage conditions must be met:

- Temperature range: -10 – 55 °C
- Max. humidity: max. 80 % (use desiccant)

9 Disposal

At the end of its service life, the UV lamp for the InnuPure C16 *touch* and all its electronic components must be disposed of in accordance with the applicable regulations as electronic waste.

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© Copyright 2017, Analytik Jena AG